METHODSNowTM
A CAS SOLUTION

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CAS[®]
A DIVISION OF THE
AMERICAN CHEMICAL SOCIETY

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
<https://www.methodsnow.com/>

METHODSNOW™
A CAS SOLUTION

Your single source to search and compare the latest published scientific methods.
Log In with your CAS Username

Learn more

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Username
methodsnow001

Password
.....

Log In

Forgot Username or Password?

Use of this product means that you agree to the Terms & Conditions.

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Analysis of Flavonoids in Trichosanthes cucumerina by Solvent extraction

CAS MN: 1-131-CAS-55818

Method Category: Natural Product Isolation Analysis
Technique: Spectrophotometry; Colorimetry; Solvent extraction

Materials	Role
Flavonoids	analyte
Leaf	matrix
Stem	matrix
Fruits	matrix
Trichosanthes cucumerina	matrix
Whatman no.1 filter paper	material
Methanol	reagent
Potassium acetate	reagent
Aluminum chloride	reagent

Materials

CAS No. 7446-70-0
AlCl₃
Aluminum chloride

Cl[Al](Cl)Cl

Close

[View Structure](#)
67-56-1

[View Structure](#)
127-08-2

[View Structure](#)
7446-70-0

Source

Total phenolics, flavonoids and antioxidant activity of *Trichosanthes cucumerina* Linn
Choudhary, Soniya; Tanwer, Babeet Singh; Vijayvergia, Rekha
Drug Invention Today (2012), 4 (5), 368 - 370. Drug Invention Today
CODEN: DIITRCZ | ISSN: 09757619

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Abstract ~
Trichosanthes cucumerina Linn. (Snake gourd) of Family: Cucurbitaceae was analyzed for their antioxidant activity as well as their total phenolic content and flavonoid contents using commonly accepted methods. The plant parts (stem, leaf and fruit) were extracted in methanol. The total level of phenolic contents (32.2 ± 0.49 mg GAE/gm DW) and flavonoids (7.82 ± 0.67 mg QE/gm DW) were found higher in leaves than other plant parts. The antioxidant activity was measured by DPPH radical scavenging activity and maximum activity was found in leaves (90.17 ± 0.67 %) in 100µg concentration and it is dose dependent.

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Equipment Used
Spectrophotometer

Conditions
Instrument
Wavelength: 415 nm

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Solvent extraction

1. Collect the stem, leaves and fruits of plant *Trichosanthes cucumerina*. L.
2. Dry the samples at room temperature, crush in grinder.
3. Extract the powder with methanol for 48 h.
4. Filter the extract through Whatman no.1 filter paper and appropriately dilute with methanol.

Determination of total flavonoid contents

1. Mix the plant extracts (0.5 mL) with 1.5 mL of methanol, 0.1 mL of 10% aluminum chloride, 0.1 mL of 1 M potassium acetate and 2.8 mL of distilled water and keep at room temperature for 30 min.
2. Measure the absorbance of the reaction mixture at 415 nm.
3. Prepare the calibration curve using 12.5 to 100 µg/mL of Quercetin in methanol.
4. Express the results as amount of flavonoid content (Quercetin equivalent, QE) per g dry weight.

Validation

Concentration	4.75 ± 1.14 mg QE/g dry weight (Stem)
	7.82 ± 0.67 mg QE/g dry weight (Leaf)
	2.05 ± 0.49 mg QE/g dry weight (Fruit)

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10

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Compare (3/3)

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	1	2	3
Title	Analysis of Flavonoids in <i>Trichosanthes cucumerina</i> by Solvent extraction	Analysis of Flavonoids in <i>Salacia chinensis</i> by Solvent extraction	Analysis of Flavonoids in <i>Prunus persica</i> by Fractionation
CAS Method Number	1-131-CAS-55818	1-131-CAS-119913	1-131-CAS-140982
Method Category	Natural Product Isolation Analysis	Natural Product Isolation Analysis	Natural Product Isolation Analysis
Technique	Spectrophotometry; Colorimetry; Solvent extraction	Spectrophotometry; Solvent extraction	Fractionation; Colorimetry; Solvent extraction
Analyte	Flavonoids	Flavonoids	Flavonoids
Matrix	Leaf; Stem; Fruits; <i>Trichosanthes cucumerina</i>	Fibrous materials; Fruits; <i>Salacia chinensis</i>	Fruits; <i>Prunus persica</i>
Other Materials	Whatman no. 1 filter paper	Acetone; Potassium acetate; Aluminum chloride; Whatman No. 1 filter paper	Acetone; Sodium hydroxide; 1-Butanol; Aluminum trichloride hexahydrate; Sodium nitrite; Buchner
Equipment Used	Spectrophotometer	UV-spectrophotometer, 190 double beam, Shimadzu, Japan	Spectrophotometer, Pharmaspec UV-1700, Shimadzu, Kyoto, Japan
Conditions	Instrument: Wavelength: 415 nm	Instrument: Wavelength: 415 nm	Instrument: detection wavelength: 510 nm

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Method	1	2	3
Method	Solvent extraction	Solvent extraction using acetone	Extraction and fractionation with butanol
Concentration	4.75 ± 1.14 mg QE/g dry weight (Stem), 7.82 ± 0.67 mg QE/g dry weight (Leaf), 2.05 ± 0.49 mg QE/g dry weight (Fruit)	0.20 ± 0.28 mg of rutin equivalents (RE) per gm fresh weight	28.26 ± 2.36 mg quercetin / g

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